

# Supplementary protocol for: *Design of highly functional genome editors by modeling the universe of CRISPR-Cas sequences*

## HEK293T transfection protocol

### Preparation of plasmid DNA for transfection:

#### **1. Add spacer of interest immediately upstream of the sgRNA scaffold.**

- a. Flanking regions for the spacer sequence should be:
  - i. Upstream: TGTGGAAAGGACGAAACACC
  - ii. Downstream (OpenCRISPR-1 sgRNA, 16 nt stem):  
GTTT TAGAGCTGTGTTGAAA
  - iii. Downstream (OpenCRISPR-1 sgRNA, 12 nt stem):  
GTTT TAGAGCTGGAAACAGC
  - iv. Downstream (SpCas9 sgRNA, 12nt stem): GTTT TAGAGCTAGAAATAGC
- b. KLD Enzyme Mix (NEB) is an efficient cloning approach for spacer insertion, but NEBuilder® HiFi DNA Assembly Master Mix (NEB) will also work with the appropriate input fragments.
- c. Platinum™ SuperFi II PCR Master Mix (Thermo Fisher) is a high fidelity polymerase able to generate very long amplicons, which is very useful for the whole plasmid amplifications needed for the KLD cloning method.
  - i. Reverse primer for KLD assembly: GGTGTTTCGTCCTTCCACAAG
  - ii. Forward primer for KLD assembly (OpenCRISPR-1 sgRNA, 16 nt stem):  
[spacer]GTTT TAGAGCTGTGTTGAAAAACAC
  - iii. To generate PCR products for KLD cloning use 10 ng of template plasmid in a 25 µL PCR reaction. Only 1 µL of this unpurified PCR product will be needed for the KLD reaction.
- d. NEB® Turbo Competent *E. coli* (High Efficiency) work well for cloning and plasmid propagation, but many other types of *E. coli* competent cells can also be used.

#### **2. Plasmid DNA purification.**

- a. From a fresh colony or glycerol stock of cells transformed with the nuclease and guide plasmids inoculate 500 µL of terrific broth (Difco™, BD Life Sciences) in a 2 mL deep well plate (Greiner Bio-One, 780271).
- b. Incubate for 12 - 24 hours at 37°C, shaking at 1000 rpm.
- c. Purify plasmid DNA using an endotoxin-free plasmid purification kit.
  - i. Mag-Bind® Ultra-Pure Plasmid DNA 96 Kit (Omega Bio-tek) is an option that scales well to large numbers of plasmids.
- d. Measure the concentration of the eluted plasmid DNA.

- i. QuantiFluor® dsDNA System (Promega) is an option that scales well to large numbers of plasmids.
- e. Normalize the DNA concentrations to 25 ng/μL by diluting in nuclease-free water.
  - i. Diluted DNA solutions should be stored at -20°C before use in transfection.

### Transfection of HEK293T cells:

#### **3. Thaw and expand HEK293T cells (ATCC, CRL-3216).**

- a. Standard media conditions are high glucose DMEM with 4 mM L-glutamine, 1 mM sodium pyruvate and phenol red pH indicator (Gibco), supplemented to a final concentration of 10% fetal bovine serum (Gibco) and 1X penicillin-streptomycin (Gibco).
- b. Standard growth conditions are 37°C and 5% (v/v) CO<sub>2</sub> in a humidified incubator.
- c. Cells should be passaged once after freezing before use in a transfection.

#### **4. Deliver nuclease and guide plasmids by transfection.**

- a. 24 hours prior to transfection seed cells at a density of 1x10<sup>3</sup> cells/well in 200 μL of supplemented DMEM in 96-well tissue culture-treated plates (Nunc™ Edge™, Thermo Fisher).
- b. Premix nuclease and guide plasmids.
  - i. Each transfection well should receive 100 ng of plasmid DNA: 50 ng of the nuclease plasmid and 50 ng of the guide plasmid.
  - ii. Premix 4 μL each of the 25 ng/μL stocks of nuclease and guide plasmid. (This prepares a double volume stock of the plasmid mix to allow for dead volume during the transfection.)
- c. Prepare and deliver transfection mixes.
  - i. Mix 4 μL of the nuclease and guide plasmid solution with 5 μL of Opti-MEM (Gibco).
  - ii. Mix 0.2 μL of TransIT®-2020 transfection reagent (Mirus Bio) with 4 μL of Opti-MEM (Gibco).
  - iii. Combine the DNA and TransIT®-2020 mixes.
  - iv. Incubate at room temperature for 15 - 30 minutes.
  - v. Add the entire mixture dropwise to the cells seeded 24 hours previously.
  - vi. Gently rock the plate to distribute the liposome complexes.
  - vii. Incubate the plate for 72 hours prior to harvesting for sequencing analysis.
    - 1. If the edited cells will be used for other downstream assays that require cell propagation, the cells should be passaged into larger culture vessels between 24 - 36 hours post-transfection.

## Sample preparation for NGS sequencing analysis:

### **5. Harvest cell lysates.**

- a. Prepare lysis buffer: 100 mM Tris-HCl, pH 7.5 (VWR); 0.05% SDS (VWR); 25 µg/mL Proteinase K (NEB).
- b. Carefully aspirate the culture media from each transfection well to be harvested. Then wash each well with 100 µL 1X PBS (Gibco), taking care to not disturb the cell layer.
- c. Add 25 µL of lysis buffer per well.
- d. Incubate at 37°C for 1 hour.
- e. Add 25 µL of nuclease-free water (Qiagen) per well.
- f. Harvest lysates into a 96-well PCR plate (VWR).
- g. Collect lysates to the bottom of the well by centrifugation (1000 xg, 2 minutes).
- h. Boil lysates at 98°C for 15 minutes in a thermocycler.

### **6. Prepare and pool NGS amplicons.**

- a. Use 5 µL of cell lysate (approximately  $1 \times 10^3$  harvested cells) as template for the target specific PCR reaction (25 µL total volume, 25 amplification cycles, 60°C primer annealing) (Invitrogen™ Platinum™ SuperFi II PCR Master Mix, Thermo Fisher), which will amplify the region of interest from the genome and append adapter sequences to allow for sample barcoding.
- b. Dilute the PCR products 1:100 in nuclease-free water (Qiagen).
- c. Use 8 µL of the diluted PCR product as template for the barcoding PCR (25 µL total volume, 12 amplification cycles, 60°C primer annealing) (Q5® High-Fidelity DNA Polymerase, NEB) which will append barcode sequences and Illumina-compatible adapters.
  - i. Barcoding primers: xGen UDI 10nt Primer Plates 1-16 (IDT).
- d. Pool 10 µL of each barcoded PCR product and mix thoroughly.
- e. Submit pooled library for NGS analysis.